# Levenhuk 1ST, Levenhuk 2ST, Levenhuk 3ST Stereo Microscopes


User Manual  
Návod k použití  
Bedienungsanleitung  
Guía del usuario  
Instrukcja obsługi  
Инструкция по эксплуатации

*Radost zaostřit  
Zoom ran und hab Fun!  
Amplie y disfrute  
Radość przybliżania  
Приближает с удовольствием*

**levenhuk**<sup>o</sup>  
Zoom&Joy


Levenhuk 1ST

EN

- 1 Pole
- 2 Locking knob
- 3 Focus knob
- 4 Eyepieces
- 5 Diopter adjustment ring
- 6 Binocular head
- 7 Objective(s)
- 8 Stage clips
- 9 Stage
- 10 Reflected light illumination
- 11 Battery compartment and ON/OFF switch
- 12 Illumination and power switch
- 13 Frosted glass stage inset
- 14 Transmitted light illumination

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
Please read this manual carefully before using your microscope for the first time.

Levenhuk 1ST, Levenhuk 2ST and Levenhuk 3ST microscopes were designed to be used during assembly of intricate equipment, including electronics; for school experiments and scientific research. All models may be found in schools and research institutes, at home and in manufacturing plants.

## Specifications

	Levenhuk 1ST	Levenhuk 2ST	Levenhuk 3ST
Head	binocular, vertical	binocular, vertical	binocular, inclined at 45°
Optics material		glass	
Magnification, x	20	40	20–40
Eyepiece tube diameter, mm	23.2	30.5	30.5
Eyepieces		WF10x (2pcs)	
Objectives	2x	4x	2x, 4x
Working distance, mm	65	60	57
Interpupillary distance, mm	55–75	60	60
Stage	50mm, with stage clips, white stage plate	with stage clips, double-sided (black and white) stage plate	matte glass stage with stage clips, additional double-sided (black and white) stage plate
Eyepiece diopter adjustment	+	+5/-5D	+5/-5D
Focus	coarse, 35mm	coarse, 40mm	coarse, 40mm
Power supply	2 AA batteries (not included in the kit)	-	220V/50Hz or 110V/60Hz (optional)
Illumination	upper LED	natural illumination	upper and lower halogen
Operating temperature range		-5... 40 °C (23... 104 °F)	
Additional accessories	dust cover	rubber eyecups	rubber eyecups; backup halogen lamp 12V/10W

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

## Using the microscope

### Illumination (1ST and 3ST Levenhuk Microscopes)

Plug your microscope into a power source or insert batteries, turn it on and choose the required illumination source. If your microscope is fitted with brightness control, you can also choose the brightness level of illumination.

### Stage inset glass

- Frosted glass stage inset may be used for observations of transparent specimens or thin slices of objects. Place it on the stage and secure it in place with the locking screw. Lower illumination source should be used during observations with this glass stage inset.
- The kit of your new microscope may also include a black and white glass stage inset. Normally, you would want to use the white side; however, if a specimen is very bright, or colored white, you should use the black side for additional contrast. Upper illumination source should be used during observations with this glass stage inset. Before installing the black and white glass make sure to remove the frosted glass from the stage.

## Preparations

Place a specimen in the center of a blank slide, put the slide on the stage and lock it in place with stage clips.

## Eyecups

For additional comfort during observations you can use rubber eyecups on both eyepieces.

## Focusing, diopter adjustment, interpupillary distance

Place a slide on the stage and lock it in place with stage clips. Loosen the locking knob on the microscope pole and slide the head up or down into a desired position and retighten the locking knob. Look through the right eyepiece and turn the objective lens handle until you see the specimen in your field of view. Afterwards, turn the focus knob until the image is sharp and clear. Look through the left eyepiece and rotate the diopter adjustment ring until the image is sharp and clear. Do not adjust focus during diopter adjustment! Upon completing this process, look through both eyepieces, hold them with your hands and bring them closer or move them apart until a double image turns into a single image.

## Changing objective lenses

Levenhuk 3ST has a revolving nosepiece with two objective lenses. To switch between them you have to rotate the nosepiece.

## Changing lamps and fuses

**CAUTION:** Before changing a lamp or a fuse make sure the microscope is unplugged from a power source. Also, wait until the lamp cools down.

## Changing the upper illumination lamp

Loosen the locking screw on the upper illumination housing and remove the light socket from within. Replace the old lamp with a new one and put the light socket back. Retighten the locking screw.

## Changing the lower illumination lamp

Loosen the locking screw on the stage and remove the glass. Remove the old lamp from the light socket and replace it with a new one. Put the glass back in its place and retighten the locking screw.

## Changing the fuse

Fuse box is located on the bottom of the microscope base. Remove the cover and replace the old fuse with a new one.

## Care and maintenance

- **Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.**
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.

- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.
- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type.
- Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit.
- **Seek medical advice immediately if a small part or a battery is swallowed.**

## Battery safety instructions

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Always replace the whole set of batteries at one time; taking care not to mix old and new ones, or batteries of different types.
- Clean the battery contacts and also those of the device prior to battery installation.
- Make sure the batteries are installed correctly with regard to polarity (+ and -).
- Remove batteries from equipment that is not to be used for an extended period of time.
- Remove used batteries promptly.
- Never attempt to recharge primary batteries as this may cause leakage, fire, or explosion.
- Never short-circuit batteries as this may lead to high temperatures, leakage, or explosion.
- Never heat batteries in order to revive them.
- Do not disassemble batteries.
- Remember to switch off devices after use.
- Keep batteries out of the reach of children, to avoid risk of ingestion, suffocation, or poisoning.
- Utilize used batteries as prescribed by your country laws.

## International Warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship. Lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase. The warranty entitles you to free repair or replacement of the Levenhuk product in any country where a Levenhuk office is located if all warranty conditions are met. For further details please visit our web site: <https://www.levenhuk.com/warranty>.

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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